
# HAND AND FOOT SURGERY

## **SNAP OFF SCREW SYSTEM**

Self-drilling and Self-tapping

## **HBS COMPRESSION SCREW SYSTEM**

Double threading for optimal compression

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### Introduction

For the fixation of intra-articular fractures the new HBS system offers a choice of Standard or High Compression Screws.

Being cannulated, the screws can be inserted over a 1 mm Guide Wire, thus eliminating the need to use a Jig, and allowing for percutaneous insertion. The self-tapping screw requires only a single cannulated drill, and the T-Drive system ensures complete control.

Since the screws are both inert and non-protrusive, they do not have to be removed, making them the ideal implant for use within or adjacent to a joint.

### Indications

- Scaphoid Fractures
- Carpal Fractures and Nonunions
- Metacarpal Fractures
- Distal Radial Fractures (articular fragments)
- Ulnar Styloid Fractures
- Radial Head Fractures
- Capitellum Fractures
- Humeral Head Fractures
- Glenoid Fractures
- Inter-Carpal Fusions
- Inter-Phalangeal Fusions
- Metatarsal Osteotomies
- Tarsal Fusions
- Malleolar Fractures
- Patellar Fractures
- Osteochondral Fractures
- Odontoid Fractures
- Mandibular Fractures

### Advantages

- Cannulated Screw for 1 mm Guide Wire
- Two kind of compressions
- Self retaining T-Drive
- Both threads are self-tapping
- Screw length in 1 mm increments

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# INSTRUMENTS

Cat. No. 91-030.01  
Graphic Case (Tray in Steel) for Set Cat. No. 91-030.00

Cat. No. 91-030.80 Sterilization Container 310x190x65 mm,  
Lid and Bottom perforated

Recommended Sterilization Container  
for Set Cat. No. 91-030.00 (are not include in the Set)

## INSTRUMENT / IMPLANT SET

Cat. No. 91-030.00 HBS and Snap off Screw  
Implant and Instrument Set

## IMPLANTS

For Cat. No. see page 6!

Scaphoid bone screws in Titanium DIN ISO 5832-3  
HBS bone screw system

Dia. 4.0 mm  Dia. 3.0 mm

For Cat. No. see page 9!

Self-drilling and tapping Snap off Screw Dia. 2.0 mm  
in Titanium DIN ISO 5832-3 (weil osteotomy)

 Dia. 2.0 mm

Listing for Set Cat. No. 91-030.00

Cat. No.: Instruments for HBS – Snap off Screw System	Pcs.:
91-030.39 HBS – Guide Wire Dia. 1.0 mm x Length 80 mm	-4-
91-030.40 HBS – Cannulated Drill Bit Dia. 2.1 mm / 3.3 mm	-1-
91-030.41 HBS – Cannulated Drill Bit Dia. 2.1 mm	-1-
91-030.42 HBS – Screw Length Gauge	-1-
91-030.43 HBS – Measuring Sleeve for Guide Wire	-1-
91-030.44 HBS – Cannulated Screw Driver, Hexagonal 2.0 mm	-1-
91-030.45 HBS – Screw Driver for Snap off Screws	-1-
95-193.00 Screw Forceps	-1-
Cat. No.: Dia. 3.0 mm Ti. HBS Screws, Cannulated	Pcs.:
91-030.12 Length 12 mm, Titanium	-4-
91-030.14 Length 14 mm, Titanium	-4-
91-030.16 Length 16 mm, Titanium	-4-
91-030.18 Length 18 mm, Titanium	-4-
91-030.20 Length 20 mm, Titanium	-4-
91-030.22 Length 22 mm, Titanium	-4-
91-030.24 Length 24 mm, Titanium	-4-
91-030.26 Length 26 mm, Titanium	-4-
91-030.28 Length 28 mm, Titanium	-4-
91-030.30 Length 30 mm, Titanium	-4-
Cat. No.: Dia. 2.0 mm Ti. Snap off Screws	Pcs.:
91-030.60 Length 11 mm, Titanium	-4-
91-030.62 Length 12 mm, Titanium	-4-
91-030.64 Length 13 mm, Titanium	-4-
91-030.66 Length 14 mm, Titanium	-4-

# INSTRUMENTS

Cat. No. 91-030.39 HBS-Guide Wire Dia. 1.0 mm x Length 80 mm

Cat. No. 91-030.40 HBS-Cannulated Drill Bit Dia. 2.1 mm/3.3 mm

Cat. No. 91-030.41 HBS-Cannulated Drill Bit Dia. 2.1 mm

Cat. No. 91-030.42 HBS-Screw Length Gauge

Cat. No. 91-030.43 HBS-Measuring Sleeve for Guide Wire 91-030.39

Cat. No. 91-030.44 HBS-Cannulated Screw Driver, Hex. 2.0 mm

Cat. No. 91-030.45 HBS-Screw Driver for Snap off Screws

Cat. No. 95-193.00 Screw Forceps

Cat. No. 95-290.00 Bone Clamp

Cat. No. 06-121.03 Inge 16.0 cm


# DOUBLE THREADING FOR OPTIMAL COMPRESSION

## Design rationale & Main features

The Compression Screw is easy to insert (over a guide wire) and provides efficient compression (through two separate threadings with different pitches, and an intermediate unthreaded section), thus ensuring quick, dependable internal fixation.

### INDICATIONS



- Distal & proximal metatarsal osteotomies
- SCARF osteotomy
- Uni and biocortical internal fixation (ex.: scaphoid)
- Small bone fusion

Cat. No.:	Total Length	Cat. No.:	Total Length
	Cannulated		non Cannulated
91-030.12	12 mm	91-030.13	12 mm
91-030.14	14 mm	91-030.15	14 mm
91-030.16	16 mm	91-030.17	16 mm
91-030.18	18 mm	91-030.19	18 mm
91-030.20	20 mm	91-030.21	20 mm
91-030.22	22 mm	91-030.23	22 mm
91-030.24	24 mm	91-030.25	24 mm
91-030.26	26 mm	91-030.27	26 mm
91-030.28	28 mm	91-030.29	28 mm
91-030.30	30 mm	91-030.31	30 mm


# SURGICAL TECHNIQUE (SCARF OSTEOTOMY)

## EXPOSURE AND EXOSTOSECTOMY

- After lateral freeing of the base of the phalanx, a medial skin incision is made over the first metatarsal.
- Exostosectomy is performed using an oscillating saw, taking care to preserve cartilage integrity.
- Edges of the cut are smoothed off using a reamer or a small rasp.

## OSTEOTOMY

- The longitudinal cut is performed on the medial aspect of the metatarsal shaft, parallel to the plantar surface.
- Transverse bone cuts should be parallel to each other, and between 45° and 60° (depending on the technique used) to the longitudinal bone cut.

## TRANSLATION


- After translation has been performed, it is maintained with the special bone clamp.
- Lateral translation is normally used. However, certain corrections may require translation in the frontal or sagittal plane (for lowering or shortening).

## FIXATION


- a. A 10/10 Kirschner wire is inserted at the proper entry point and with the proper angulation (for head or shaft fixation), to serve as a guide for later drilling and screw insertion.
- b. Use the screw length gauge (using the subtraction principle) to determine the appropriate length of the screw. *The lag screw should be at least 4 mm shorter than the measured length to avoid cartilage penetration.*
- c. The cannulated drill is inserted over the guide wire and fully advanced to create the countersink for the screw head.
- d. The selected screw is inserted and its head is carefully countersunk to ensure optimal compression and avoid later impingement. *Make sure that the diaphyseal screw is firmly anchored in both cortices.*

The proximal screw is inserted using the same technique.

## ANTEROMEDIAL RESECTION

- Once the screws are positioned, the anteromedial angle is cut in line with the exostosectomy, using an oscillating saw. Edges of the cut are smoothed off.
- The capsule is closed in a routine fashion.

# SELF-DRILLING & SELF-TAPPING

## Design rationale & Main features


The Snap off Screw provides superior fixation: it saves time (no need for a pilot drill hole), and it is easy to use, safe (clean break), accurate (guide point), and efficient (self-compression).

The Snap off Screw consists of two parts: implantable screw which provides firm anchorage inserter which allows powered insertion.

### INDICATIONS

- Weil osteotomy
- Unicortical internal fixation

- 91-030.60 Snap off Screw Ø 2 mm Length 11 mm
- 91-030.62 Snap off Screw Ø 2 mm Length 12 mm
- 91-030.64 Snap off Screw Ø 2 mm Length 13 mm
- 91-030.66 Snap off Screw Ø 2 mm Length 14 mm

## SURGICAL TECHNIQUE (WEIL OSTEOTOMY)

### EXPOSURE

The procedure is performed using a dorsal intermetatarsal and/or transverse approach. After the two extensor muscles have been separated:


- Hohmann retractors are placed on both metatarsal sides.
- The metatarsophalangeal joint is dislocated between the extensor digitorum longus and the extensor digitorum brevis.
- A Hinge spreader is inserted to protect the extensor muscles and afford good exposure for the osteotomy.

## OSTEOTOMY

Osteotomy is performed using an oscillating saw:

- Make a 3 cm (approximately) horizontal cut parallel to the sole, to increase the interfragmental contact area and thus enhance healing.
- Osteotomy results in spontaneous recession of the metatarsal head, which relieves tension on soft tissue.


## TRANSLATION

- Grasp the metatarsal head with Kocher forceps.
- Use the „Index Plus Minus“ formula and the Lelièvre Curve to determine the amount of recession of the metatarsal head.
- The metatarsal head must be held in the correct position for subsequent screw fixation.


## INSERTION OF THE SNAP OFF SCREW

- Connect the screw inserter to the power drill, and drive the screw into the metatarsal.
- The inserter snaps off as soon as the screw head makes contact with the dorsal cortex.
- If necessary, insertion of the screw can be completed with the special screwdriver (with 3 notches).


## RESECTION OF THE BONE PEAK

- Bone peak is resected using Liston pliers. This allows deep flexion of the metatarsophalangeal joint.
- It may be necessary to perform a Z-shaped release (Green technique) of the extensor muscles.

## NOTICE

## Orthomedical GmbH Implants

Hauptstraße 5  
D - 78589 Dürbheim  
Tel. +49 (0) 7424/9403 - 24  
Fax +49 (0) 7424/9403 - 43

[www.orthomedical.de](http://www.orthomedical.de)  
[mail@orthomedical.de](mailto:mail@orthomedical.de)