**Oszillierende Gipssäge**  
**Oscillo-Saw**  
**Scie oscillatrice**  
**Sierra Oscillo**  
**Sega oscillante**

Gebrauchsanweisung / Operating Manual  
Notice d'utilisation / Instrucciones de  
servicio / Istruzioni per l'uso

## Oszillosäge / Oscillo-Saw / Scie oscillatrice / Sierra Oscillo/ Sega oscillante

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Seite 3

DEUTSCH

### **Operating Manual**

Page 15

ENGLISH

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Page 27

FRANÇAIS

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Página 39

ESPAÑOL

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Pagina 51

ITALIANO

DEUTSCH

# Oszillierende Gipssäge

## Inhalt

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# Oscillating cast cutting saw

## Content

1	Scope of delivery .....	17
2	Safety remarks .....	17
3	Mounting / exchanging the saw blade .....	18
4	Operation .....	19
5	Cleaning and maintenance .....	21
6	Spare Parts .....	21
7	Servicing .....	21
8	Warranty .....	25
9	Environmental protection .....	25
10	Specifications .....	25
11	Manufacturer and aftersales service address .....	25

## Oscillating cast cutting saw

## 1 Scope of delivery

Oscillating cast cutting saw
Circular saw blade, dia. 50 mm, x-bite
Circular saw blade, dia. 65 mm, x-bite
2x fork wrenches
Operating instructions

## 2 Safety remarks

All safety remarks in these operating instructions are marked with the symbol **!**.

**Before using the saw, carefully read through these operating instructions and the valid national occupational safety regulations and act accordingly.**  
*Please keep these operating instructions for later reference and always include them if the device is sold or passed on to third parties.*

**Application in accordance with the intended application:  
This device is intended for sawing open plaster casts.**

**Never rivet or screw signs or symbols onto the device, as this can render the protective insulation ineffective. We recommend using adhesive signs.**

**Only ever use undamaged plugs and cables.  
Check cables and plugs regularly.**

**The mains voltage and specified voltage rating shown at the device must be in agreement.**

**Only ever use original accessories.**

**Wear personal protective gear such as protective goggles,  
ear protectors and gloves.**

**The device must not be operated in areas at risk from explosion.**

### 3 Mounting / exchanging the saw blade

**Danger of injury due to unintentional activation.  
Pull out the plug before changing the saw blade.  
If applicable use thick gloves to prevent injury when handling the saw blade.  
Only ever connect the oscillating cast cutting saw to the mains when it is  
switched off.**

#### 3.1 Releasing the retaining screw / removing the old saw blade

1. Place the cast cutting saw on a stable underlay (table).
2. Apply one fork wrench at the drive shaft behind the saw blade, and the other one at the front at the retaining screw (see fig. 1).
3. Hold the rear wrench firmly and use it to fix the cast cutting saw.
4. Release the retaining screw using the front wrench.
5. Remove the retaining screw and the fixing ring. Ensure that you do not lose any components.
6. If applicable, remove the old saw blade.

**Fig. 1:**  
Releasing the retaining screw

### 3.2 Inserting the new saw blade

1. Check that the female thread and saw blade support of the drive shaft are clean.
2. Position the new saw blade. *It makes no difference which side of the saw blade is facing inwards or outwards.*

**Important!**

The notches in the saw blade must be positioned precisely flush on the spigots of the drive shaft.

3. Position the fixing ring. The notches in the fixing ring must always be pointing towards the saw (see fig. 2). These must also be positioned precisely on the spigots.
4. Screw the retaining screw in manually until hand tight.
5. Then tighten firmly using the fork wrench, as described above. The tightening torque is 6 - 7 Nm. Caution! Do not overtighten.

**Fig. 2:**  
Inserting the saw blade

## 4 Operation

### Switching on and off (see fig. 3)

- Switching on: Push the switch forwards ( I )  
Switching off: Push the switch back ( 0 )

Fig. 3

### Setting the oscillating frequency (see fig. 4)

Select the oscillating frequency when the motor is running. The optimum oscillating frequency can be selected at the speed gauge depending on the properties of the plaster cast.

Oscillating frequency preselection using the electronic speed gauge.

- Lowest oscillating frequency ( 1 )  
Highest oscillating frequency ( 8 )

Fig. 4

## 4.1 Opening plaster casts

First select a suitable saw blade.

Apply the oscillating saw with the saw blade mounted to the cast you wish to open and press in lightly. Your guiding hand can be used as a support and for depth control, and to prevent sudden penetration of the blade when the cast has been cut through.

As long as the saw blade is cutting through the plaster, you will be able to feel a resistance. Once the plaster cast has been cut through, this resistance disappears. Once the breakthrough point is tangible, lift the saw blade off lightly without removing it from the cutting groove and move it forwards appr. 15 mm in the direction of cutting.

Then apply the blade again and proceed as described above until the next piece has been cut through.

In this way, the plaster cast is more safely opened through a series of cuts and linear movements along the cutting line than would be the case if executing a continuous linear cut.

After you have used the saw only a few times, you will become accustomed to guiding the saw through the plaster cast and you will also develop an instinctive feel for the moment where no further resistance is tangible and no further pressure should be applied to the saw blade.

If the saw blade oscillates for too long on the wadding, high-speed oscillation can bring about a burning sensation on the skin under the wadding.

### Additional remarks:

- Because different materials are used, every plaster or plastic cast has different properties. While working, try to determine the optimum speed for the respective material.
- A lower speed is advantages when cutting through plastic casts. This prevents the plastic from starting to melt due to friction heat generated during cutting.
- Where possible, avoid cutting through plaster casts without wadding when using the oscillating cast cutting saw. Injuries can result if the skin has become stuck to the plaster, and so does not oscillate with the saw blade. Where a plaster cast close to the skin is preferred, we recommend laying a knitted or felt strip along one side of the arm or leg and drawing a line precisely above this strip on the surface of the plaster cast using indelible ink. When sawing open the cast, this line can be used as guiding mark.
- When applying a plaster cast following surgery, the position of the operation wound should be marked on the surface of the plaster cast. This will allow a window to be cut subsequently in the plaster for removing stitches and observation of the wound.
- Cutting a window in the plaster also simplifies the extraction of bone nails and also facilitates any necessary punctures or incisions.
- Metal rails inserted in the plaster cast should not be cut.
- While serving, check the retaining screw for a firm fit and if necessary tighten.
- If the motor is overloaded, the device is automatically switched off. The fuse (see fig. 9 on page 24, item 110) is located on the electronic circuit board (see fig. 7 on page 24, item 110) and can be quickly and easily exchanged after unscrewing the housing.
- Used and blunt saw blades place excessive stress on the motor and should consequently be continuously checked and exchanged in good time.
- To ensure even wear of the saw blade, open the retaining screw and turn the saw blade by 90°.

## 5 Cleaning and maintenance

**Danger of injury through unintentional activation.  
Pull out the plug before any cleaning or main-tenance work.**

The following measures should be performed once a week or more often if the saw is used frequently: Clean the cooling slots (see fig. 5) to ensure sufficient ventilation for cooling the motor. Clean the motor compartment from the outside with a jet of dry compressed air.

Fig. 5

## 6 Spare Parts

Fork wrench

Retaining screw


Fixing ring

## 7 Servicing

### 7.1 Saw head

(Drawing items refer to fig. 6)

Pos.	Pc.	Description
10	1	Saw head housing
20	1	main shaft
30	1	Floating bearing front
40	2	Head Screw with washer
50	1	Floating bearing back
70	1	Fork
80	1	Fork screw
90	2	Bearing screw


**Fig. 6:**  
Explosion view of the saw head


## 7.2 Saw motor

The drawing items refer to fig. 7 to 9, next page.


Pos.	Pc.	Description
10	1	Motor housing
20	1	Switching rod
30	1	Switching slide
40	1	Pole housing with winding
50	2	Ejot Delta PT screw
60	1	Air guiding ring
70	2	Brush holder
75	2	Plug connector
80	2	Carbon brush
90	1	Contact spring
110	1	Electronic circuit
115	1	Fuse
116	1	Pressure pad for fuse
130	1	Cable clamping bridge
150	1	Lid
155	1	Pressure pad for lid
160	1	Ejot Delta PT screw
170	5	Ejot Delta PT screw
190	1	Protective hose
250	1	Anchor
260	1	Deep groove ball bearing
270	1	Deep groove ball bearing
300	1	Intermediate bearing
310	1	Bearing bush
320	1	Deep groove ball bearing
330	4	Ejot Delta PT screw
340	1	Parallel pin
180	1	Cable standard "Europe", 3-pin
180	1	Cable "USA", 2-pin
180	1	Cable "Australia", 3-pin
180	1	Cable "UK", 3-pin
180	1	Cable "UK", 3-pin
180	1	Cable "Switzerland", 2-pin

**Fig 7:**  
Motor drawing

**Fig 8:**  
Motor terminal plan

**Fig 9:**  
Electronic circuit

## 8 Warranty

We provide a guarantee the cast cutting saws in accordance with statutory and country-specific legislation (proof of purchase by invoice or delivery note). The minimum warranty period is 18 month. In the event of material or production defects, we provide free repair of the device. Any attempt to carry out unauthorized repairs will result in forfeiture of your warranty cover.

Damage caused by incorrect handling, overloading or normal wear and tear is not covered by the warranty. Please contact our aftersales service department.

In case of warranty claims and repairs, please enclose a copy of the invoice or delivery note with the device.

## 9 Environmental protection

Any packaging materials, disused devices and accessories must be sent for environmentally friendly recycling or disposal. Electrical devices may not be disposed of with the household waste but taken to the local disposal and recycling facility for electrical devices.

For more information on this subject, consult your specialist dealer.

## 10 Specifications

Model	Oscillating cast cutting saw (A / G / S)	Oscillating cast cutting saw (B / V)
Mains connection	220-240 Volt, 50/60 Hz.	110-120 Volt, 50/60 Hz.
Rated current consumption	250 Watt	250 Watt
Current type	1 ~	1 ~
Oscillations	ca. 12000 - 21000 min <sup>-1</sup>	ca. 12000 - 21000 min <sup>-1</sup>
Weight	1,2 kg	1,2 kg
Protection class	II 	II 
IP Code	IP2X	
Noise level (A-graded)	Typically 75 dB (A). The noise level during sawing can exceed 85 dB (A). From 85 dB (A), wear ear protectors!	
Hand-arm vibration	Typically lower than 2.5 m/s <sup>2</sup> . Measured values determined in accordance with EN 50 144.	

## 11 Manufacturer and aftersales service address

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ORTHO-MEDICAL GmbH