

## BLACK CERAMIC COATING with SILICON COATED HANDLE INSTRUMENTS

Kerrison handle for Neuro/Spinal  
Individually made for our customers



Our CERAMIC and SILICON coated instruments offer the following advantages:

- Smooth action due to less surface friction
- Scratch resistance
- Maximum durability
- Non reflective
- Higher corrosion resistance
- Aesthetical
- Validated biocompatibility
- Easy to clean due to improved surface structure
- Universal color coding for quick identification of the bite size
- Stainless steel shafts with standard bite sizes: 1mm / 2mm / 3mm / 4mm / 5mm
- Shaft available with or without T-Slot tissue post for enhance cleaning and removal of tissue
- Special order options include: handle color, bite style, thin footplate and shaft length
- The comfortable silicone-coated handle provides improved tactile feel and increases surgeon comfort