



made in GERMANY

ORIHOMEDICAL
G m b H · I M P L A N T S




SMART MICRO SCISSORS

SMART MICRO SCISSORS


SMART MICRO FINE BLADES CUT

Available Blades



One
Serrated
Blade

Smooth
Blades



Straight

Curved

45° Angled


60° Angled

90° Angled

Angles are designed for difficult areas even far distal to molars in posterior zones.

Comfortable Handle Grip!


MICRO SCISSORS



Micro Scissor Fine Blades	
03-315.214	15cm Straight
03-315.215	15cm Straight Serrated Blade
03-315.216	18cm Straight
03-315-217	18cm Straight Serrated Blade
Micro Scissor Fine Blades	
03-315.218	15cm Curved
03-315.219	15cm Curved Serrated Blade
03-315.220	18cm Curved
03-315.221	18cm Curved Serrated Blade
Micro Scissor Fine Blades	
03-315.222	15cm 45 degree angled
03-315.223	15cm 45 degree angled Serrated Blade
03-315.224	18cm 45 degree angled
03-315.225	18cm 45 degree angled Serrated Blade

Micro Scissor Fine Blades	
03-315.630	15cm 60 degree angled
03-315.632	15cm 60 degree angled Serrated Blade
03-315.634	18cm 60 degree angled
03-315.636	18cm 60 degree angled Serrated Blade
Micro Scissor Fine Blades	
03-315.638	15cm 90 degree angled
03-315.640	15cm 90 degree angled Serrated Blade
03-315.642	18cm 90 degree angled
03-315.644	18cm 90 degree angled Serrated Blade

MICRO SCISSOR FINE BLADES 11CM




03-315.652	03.315.654	03-315.656	03-315.658	03-315.660
11cm	11cm	11cm	11cm	11cm
Straight	Curved	45 Degree	60 Degree	90 Degree

IRIS MICRO BLADE



315-622
Iris Scissor Straight 11,5cm Micro Blade



315-624
Iris Scissor Curved 11,5cm Micro Blade

