

ORLOS__AO FLEXIBLE TWIST DRILL SYSTEM

Edition: 2019/04 Revision: 03.04.2019

	Art. No.	Gesamtlänge Total Length	Nutzlänge Useable Length	Spiral-Länge Twist Length
Ø 2,5 mm	95-169.22F	35,0 mm	25,0 mm	20,0 mm
	95-169.23F	50,0 mm	40,0 mm	30,0 mm
	95-169.24F	60,0 mm	50,0 mm	40,0 mm
Ø 3,2 mm	95-169.32F	35,0 mm	25,0 mm	20,0 mm
	95-169.33F	50,0 mm	40,0 mm	30,0 mm
	95-169.34F	60,0 mm	50,0 mm	40,0 mm
	95-169.35F	70,0 mm	60,0 mm	40,0 mm
Ø 3,5 mm	95-169.42F	35,0 mm	25,0 mm	20,0 mm
	95-169.43F	50,0 mm	40,0 mm	30,0 mm
	95-169.44F	60,0 mm	50,0 mm	40,0 mm
	95-169.45F	70,0 mm	60,0 mm	40,0 mm
Ø 4,5 mm	95-169.53F	45,0 mm	35,0 mm	30,0 mm
	95-169.54F	60,0 mm	50,0 mm	40,0 mm
Bohrer auswechselbar / drills changeable				
Flexibler Schaft mit Schnellkupplung Flexible Shank, quick coupling				
	95-169.00F	135,0 mm		

