

made in GERMANY



ORLOS IN SITU PIN SHEAR CUTTER, 29 CM 30-137.29



Shear 3 different pin diameters with 1 instrument.

- Double action design allows for ease in cutting.
- Grooved handles are ergonomically designed for a comfortable and secure grip.
- Shears 2.5mm [.098"] 3.0mm [.118"] and 3.5mm [.138"] diameter pins with ease and provides a clean cut without burrs or sharp edges.



Edition 2020|01 | Revision: 02.01.2020
Printed in Malaysia | Copyright by ORTHO MEDICAL GmbH
We reserve the right to make alterations