

made in GERMANY



## ORLOS SILICONE IN SITU ROD SHEAR CUTTER, WITH DETACHABLE HANDLES, OAL=57CM 30-136.57

- Rods are sheared, leaving a smooth and clean surface rather than a sharp, jagged, burr - like surface common when using a standard pinching - type rod cutter.



Rod cut by standard pinching-type rod cutter



Rod cut by ORLOS Silicone In Situ Cutter

The ORLOS Rod Shears easily cuts: 4.5mm [.177"], 5.5mm [.217"] and 6.35mm [.250"] diameter pins and rods.



made in GERMANY



## ORLOS SILICONE IN SITU ROD SHEAR CUTTER, WITH DETACHABLE HANDLES, OAL=57CM 30-136.57

detached view



• Double action design allows for ease in cutting.

• Quick handle release

• Available with detachable handles.



• Note: Instrument is not designed to be used inside the patient's body.

Edition 2020|01 | Revision: 02.01.2020  
Printed in Malaysia | Copyright by ORTHO MEDICAL GmbH  
We reserve the right to make alterations