


## THE SPINE ENDOSCOPE



**The perfect endoscopic view**

## THE SPINE ENDOSCOPE



### Features

- Sapphire glass window
- Gold soldered opto-mechanics
- Large working channels
- Incorporated irrigation and suction channels
- Stainless Steel
- Fully autoclavable



Specification/ Item Number	30-0734-00	30-0780-00
Viewing Direction	30°	30°
Outer Diameter	Ø 6,3 mm	Ø 6,9 mm
Working length	193mm	193mm
Total length	260mm	260mm
Working channel	Ø 3,75 mm	Ø 4,1 mm
Suction/ irrigation channel	2x Ø1,4mm	2x Ø1,4mm

Cat. No.	Description
30-0734-00	Spine Endoscope Ø 6.3 mm L= 193mm, 30° Autoclavable 134°
30-0780-00	Spine Endoscope Ø 6.9 mm L= 193mm, 30° Autoclavable 134°

Edition 2021|12 | Revision: 21.12.2021  
Copyright by ORTHO MEDICAL GmbH  
We reserve the right to make alterations